



Data Distributing
the Image Distribution People



UVC SANITIZING CABINET

CLEAN™

Rapid disinfection of laptops, tablets, smartphones, call devices, noninvasive portable clinical items & more



UVC Germicidal



Safe to Operate



360° Sanitization



- Chemical-free UV-C surface disinfection process
- Effectively kill 99.9% of infectious pathogens like influenza, the common cold and coronavirus
- 360° reach - eight lamps & reflective lining
- Configurable, medical grade stainless steel racks

- Customizable disinfection cycle times
- Intuitive, easy to use and maintain
- UV absorbing, transparent window with electro-magnetic locking system
- Lifetime frame warranty



Support infection control initiatives with ultraviolet (UVC) light sanitization to encourage more frequent disinfection without the use of corrosive chemicals. Clean™ effectively kills 99.9% of infectious pathogens from laptops, tablets, smartphones, call devices, noninvasive portable clinical items & more. The UV absorbing, transparent window with electro-magnetic locking system keeps users safe while devices are being disinfected.

Designed with the highest level of flexibility and versatility in mind, Clean™ offers the ability to configure medical grade stainless steel racks to adapt to devices of varying shapes and sizes. Disinfection cycle times, from 2 minutes to 20 minutes, can be customized depending on workflow and efficacy needs.

Fitment & Dimensions

Number of Devices	Exterior Dimensions		Interior Dimensions		Product Weight	Shipping Weight
12 - 15.24"x 20.375"	H	19.7" 50 cm	H	15.25" 38.7 cm	75 lbs	90 lbs
(Maximum Device Size)	W	25" 63.5 cm	W	16.25" 41.2 cm	34 kg	40 kg
	D	21.7" 55 cm	D	20.37" 51.7 cm		

Clean™ sanitizing cabinet has been independently tested by SGS* for antimicrobial effectiveness and achieves the disinfected efficiency to 99.99% kill infectious pathogens within 5 minutes or less reducing the risk of spreading common bacteria and viruses.

SGS Test Results

Bacteria / Virus				
Escherichia coli (CICC 10389)	Staphylococcus aureus (ATCC 6538)	Pseudomonas aeruginosa (ATCC 9027)	Candida albicans (CICC 1965)	SARS-Cov-2 virus
Kill rate: 99.99%	Kill rate: 99.99%	Kill rate: 99.99%	Kill rate: 99.99%	Kill rate: 99.99%
Disinfecting time: 5 Minutes	Disinfecting time: 5 Minutes	Disinfecting time: 5 Minutes	Disinfecting time: 5 Minutes	Disinfecting time: 5 Minutes

Laboratory Information: Testing conducted by SGS's ultra-trace industrial safety and hygiene laboratory.

Models

SKU	Description
PT-CLEAN	12 Bay Configurable UV-C 360 Cabinet - 110V NA
PT-CLEAN-UK	12 Bay Configurable UV-C 360 Cabinet - 220V UK
PT-CLEAN-EU	12 Bay Configurable UV-C 360 Cabinet - 220V EU

This product is designed to conform to Electromagnetic Compatibility (EMC) standard IEC 60601-1-2 and will operate correctly in combination with other medical equipment which also meets this requirement. To avoid interference problems affecting the device, connect only equipment that complies with IEC, ISO, UL/ANSI, or CSA standards that are relevant to that equipment.

