

Process DICOM Video Files and DVD DISCs

Improve patient care by sizing, post-processing, and DICOMizing videos for archiving in PACS and VNA.

Use Vertex Video to create video snippets and single-frame images from DVD DISCs and video files. Convert it all to DICOM for PACS and VNA storage.

Process Stored DVDs

Empty the DVD cabinets. Size, edit, and convert the content for archiving in PACS and VNA.

Convert and Store Existing Video Files

Use Vertex to edit and convert popular video formats to DICOM.
3GP, ASF, WMV, AVI, MPEG, VOB(DVD)

Flexible Formatting

Store converted videos in DICOM format, compatible with your PACS/VNA (multi-frame and/or MPEG4)

Appropriately Sized Video

To minimize storage requirements Vertex Video software allows you to configure frame rate, transfer syntax, color space conversion, and auto split.

Video Processing Profiles

To minimize error and improve efficiency, multiple profiles are preconfigured for processing videos and various workflows.

Consistent DICOMization


Configure department information, and create pre-defined study and series descriptions, ensuring data integrity when creating new studies.

Client Access Control

For security and regulatory purposes, individual users, groups, or workstations can be configured to have access to incoming devices and destinations.

HIS Friendly

Add procedure codes to pre-defined study descriptions to ensure proper billing and reimbursement



Sorna